

HPS proudly serves the scientific community with high quality spectrometric standards and reference materials for the calibration of instruments, including, but not limited to, AAS, IC, ICP, ICP-MS, GC, and GC-MS.



6.0 x 6.0 inch (15.2 x 15.2 cm) 80 wipes Catalog # INT-IPR70 CW66 For wiping surfaces, fixtures and equipment

## We've Got Wipes

## LabExact<sup>®</sup> Cotton Cloth IPA Wipes

- 70% IPA (Isopropyl Alcohol) Wipes (non-sterile) meets major health and regulatory agency recommendation for concentration strength.
- Wetted sufficiently for effective use on tables, counters, carts, fixtures, door handles, and equipment.
- Excellent lay-flat characteristics. Reduced pull, tear and bunch on rough surfaces.
- Much thicker, wetter and more durable than synthetic non-woven veil materials. 100% high purity, 3x rinsed virgin cotton cloth similar to that of a handkerchief or bedsheet. 400 thread count.



- 6 x 6 inch (15.2 x 15.2 cm) dimension ideal for full hand coverage in wiping. 80 wipes per container, each easily removed with gloved hand.
- Environmentally safe cotton is 100% biodegradable and incinerates cleanly. No dyes, ammonia, bleach or fragrance added.
- Gasket seal lid prolongs product life between usages. To maintain product performance, close lid tightly after each use.
- Proudly manufactured in the USA by I.W. Tremont.
- IMPORTANT: Non-sterile, not intended for use in applications of direct patient skin contact as antiseptic or as a high level medical device disinfectant. Professional use, not intended or labeled for retail consumer sale. Batch and lot information label applied to container.

Description	Quantity	Part #	List Price
LabExact <sup>®</sup> Cotton Cloth IPA Wipes	80 wipes per container	IWT-IPA70-CW66	\$19.00

7221 Investment Drive • North Charleston, SC 29418 843.767.7900 • info@highpuritystandards.com • highpuritystandards.com